


# ALMA DENT



Catalogue 2020



With 25 years of dental experience, DENTAC serves dentists and endodontists with products designed and manufactured under **EN ISO 13485:2016**, committed to offer the experience to customers globally. Dentac offers two main lines in its portfolio as **Endodontic** and **Restorative** lines.



T-Endo must M25



T-Endo must M40



T-Endo must M50



T-Endo must tg

**T-endo MUST**, designed to provide excellent and fast root canal shaping. you can finish your treatments effectively in a short time. T-endo MUST is designed for multi-use. Thanks to the strong root canal files produced by the heat treatment technology, a short and effortless preparation in the treatment will be sufficient.

The power of T-endo MUST technology and reciprocation movement guarantees you the most effective treatment.

- Reciprocating movement to counter-clockwise
- Cutting efficiency by S-shaped cross-section
- Maximum resistance against fatigue with tm-wire technology
- Excellent material for multi-use
- One file for most of your treatments
- Packed in sterile blisters
- ISO color stoppers
- Heat treated NiTi alloy

21mm      25mm      31mm

---

T-endo MUST Reciprocal Endo File 25.06 (M25)

---

T-endo MUST Reciprocal Endo File 40.04 (M40)

---

T-endo MUST Reciprocal Endo File 50.04 (M50)

---

T-endo MUST Glide Path (TG)

---

T-endo MUST Asorted





### sealart - Resin Based Root Canal Sealer

SEALART is an epoxy resin based root canal sealer which is a paste-paste type of dual syringe. It has an outstanding chemical and physical characteristic including extremely excellent sealing property and biocompatibility

sealart Dual Syringe 13.5 gr



### bioserra - Bioceramic root canal sealer

Bioceramic sealer is used for root canal filling after the removal of infected pulp tissue in the root canal. BIOSERRA is a fluid paste and can be immediately applied inside the root canal using a disposable tip and a thickening agent. This product characteristically cures slowly by absorbing the ambient water inside the root canal. It is white and aesthetic.

bioserra Syringe 2gr

### fillart - Hydraulic temporary restorative material

Simple and fast application with filling instrument. Quick and void-free curing in a moist environment. Slight expansion of the filling material ensures a well-sealed margin

fillart 40gr



### Gutta Percha Points

Hand-rolled sterile points from highest quality Gutta Percha to extremely close tolerances. This true to size points are ideal for root canal filling after the preparation with T-endo files

### Paper Points

Highly absorbent sterile paper points for a quick drying of the root canal





Available with following shades

### triofill + - Restorative composite

Triofill+ is a **Nanohybrid** dental composite completely developed and manufactured to match both functionality and esthetics. Excellent chemical components provide improved adaptation and optimized easier handling, sculptability and polishing.

- Can be used in both anterior and posterior area
- Does not stick to composite instruments
- Half condensable, hard consistency
- Adaptable to natural tooth colors

**triofill** is also available as **Microhybrid composite**



### Transparent Matrix

Produced with transparent polyester material. It has been designed to shape the composite fillings easily. Flexible designed matrix band for ideal contact points, ensures light transmission and easy polymerisation of composite due to its transparency.

The band with 10mm width and 97mm length is compatible with all Standard matrix retainers.

- Easy to use
- Universal
- Light Transmission
- Flexible



### Articulating Paper

Special flexible paper material with a high coloring capability adapts to occlusal surfaces accurately. Thanks to its progressive structure, it enables to retain especially stressed contact points during occlusion. The paper material with only 80\_ thickness became waterproof material with hydrophilic waxes and pharmaceutical oils. One side of the paper is blue and the other side is red so that it enables to distinguish easily the stressed contact points either from lower jaw or upper jaw.

- User Friendly
- Waterproof
- Double Sided, Double Colores
- Accurate



## Prophylaxis Paste

Dentaclear Prophylaxis Paste is used for polishing of fillings and enamel after tartar removal. Mechanical abrasion is effective in enamel polishing. The pink prophylaxis material used in prophylaxis treatment leaves a pleasant taste in the mouth of the patient. The material has been manufactured with pleasant aromas that ensures patient's comfort.

- Pleasant **MINT** and **STRAWBERRY** flavor
- Easy to rinse
- Gluten free
- Economic package 340g



## UltraFace Shield

UltraFace Shield is a transparent PVC face protection against sprayed fluids, droplets, dust and other particles which can injure or even cause infections.

- Easy to wear
- Universal size
- Light weight
- Cleaning and disinfection with alcoholic solutions
- 2 UltraFace Shields in 1 box



## T-mask – surgical CLA

3-ply surgical face mask with high filter performance providing total protection and exceptional comfort. All layers consists of odourless soft material, free of glass fibres, latex and other allergenic substances. External layers water-repellent Non-Woven (PP) fabric for protection from dust and moisture. Inner layer Meltblown filter for virus and bacteria filtration.

- 3-ply with Meltblown filter
- Ultra soft
- Intergrated plastic brace
- Bacterial Filtration Efficiency (BFE) >99%
- Particle Filtration Efficiency (PFE) >99%
- Colors: green / blue / white